

### **Density Standard 2,5807g/ml @40°C (ASTM D1480)**

Art. ID	REA-DEN40270PY
Unit	100 mL
Deliverydetails	No Dangerous Good /not restricted

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 40°C	Density		2,5807	g/mL		