

Density Standard 1,1512g/ml @50°C (ASTM D1480)

Art. ID REA-DEN50170PY
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 50°C | Density | | 1,1512 | g/mL | | |