

### **Density Standard 2,0551g/ml @50°C (ASTM D1480)**

Art. ID REA-DEN50240PY  
Unit 100 mL  
Deliverydetails No Dangerous Good /not restricted

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 50°C	Density		2,0551	g/mL		