

**Density Standard 0,7572g/ml @60°C (ASTM D1480)**

Art. ID REA-DEN60080PY  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2

**Description**

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 60°C	Density		0,7572	g/mL		