

**Density Standard 0,9990g/ml @60°C (ASTM D1480)**

Art. ID REA-DEN60150PY  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1105 III E1

**Description**

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 60°C	Density		0,999	g/mL		