

## **Bingham Mortar Mixture for Rheological Measurements**

Art. ID	NIST-2493
Unit	kit for two batches
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) is intended for use in calibrating rheometers for measuring the rheological properties of cement paste and mortar. A unit of NIST-2493 consists of four containers, one glass bottle of corn syrup (500 g), two plastic jars of limestone powder (600 g each) and one plastic jar of 1 mm glass beads (1500 g) /// For certified values please see certificate