

## **Out of Stock - Item is not available at this time - C-Reactive Protein Solution**

Art. ID	NIST-2924
Unit	3 vials x 1 mL
Deliverydetails	Dry ice shipment /not restricted

### Description

This Standard Reference Material (SRM) is primarily intended for use in calibrating procedures and devices for the determination of C-reactive protein (CRP) in human serum. It can also be used for value-assignment of calibrators and control materials. A unit of NIST-2924 consists of three vials, each containing approximately 1 mL of a solution of recombinant CRP /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
k = 2,15	Certified CRP concentration		20,6 ± 1,2	µmol/L		