

ICP Multi Element Calibrat ion Standard 27 Elements each at 100 µg/mL in 2-5% HNO3

Art. ID REA-REICPCAL27DD
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1759 III

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	100	µg/mL		
	Aluminium (Al)	[7429-90-5]	100	µg/mL		
	Arsenic (As)	[7440-38-2]	100	µg/mL		
	Barium (Ba)	[7440-39-3]	100	µg/mL		
	Beryllium (Be)	[7440-41-7]	100	µg/mL		
	Bismuth (Bi)	[7440-69-9]	100	µg/mL		
	Cadmium (Cd)	[7440-43-9]	100	µg/mL		
	Cobalt (Co)	[7440-48-4]	100	µg/mL		
	Chromium (Cr)	[7440-47-3]	100	µg/mL		
	Caesium (Cs)	[7440-46-2]	100	µg/mL		
	Copper (Cu)	[7440-50-8]	100	µg/mL		
	Iron (Fe)	[7439-89-6]	100	µg/mL		
	Gallium (Ga)	[7440-55-3]	100	µg/mL		
	Mercury (Hg)	[7439-97-6]	100	µg/mL		
	Indium (In)	[7440-74-6]	100	µg/mL		
	Lithium (Li)	[7439-93-2]	100	µg/mL		
	Manganese (Mn)	[7439-96-5]	100	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	100	µg/mL		
	Nickel (Ni)	[7440-02-0]	100	µg/mL		
	Lead (Pb)	[7439-92-1]	100	µg/mL		
	Rubidium (Rb)	[7440-17-7]	100	µg/mL		
	Selenium (Se)	[7782-49-2]	100	µg/mL		
	Strontium (Sr)	[7440-24-6]	100	µg/mL		
	Titanium (Ti)	[7440-32-6]	100	µg/mL		
	Uranium (U)	[7440-61-1]	100	µg/mL		
	Vanadium (V)	[7440-62-2]	100	µg/mL		
	Zinc (Zn)	[7440-66-6]	100	µg/mL		