

Copper base, Type: Brass, RM, Chips 50 g

Art. ID DH-SL02-01
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,022	%		
	Arsenic (As)	[7440-38-2]	0,076	%		
	Bismuth (Bi)	[7440-69-9]	0,006	%		
	Copper (Cu)	[7440-50-8]	81,84	%		
	Iron (Fe)	[7439-89-6]	0,677	%		
	Manganese (Mn)	[7439-96-5]	0,035	%		
	Nickel (Ni)	[7440-02-0]	0,795	%		
	Phosphorus (P)	[7723-14-0]	0,046	%		
	Lead (Pb)	[7439-92-1]	1,17	%		
	Antimony (Sb)	[7440-36-0]	0,104	%		
	Tin (Sn)	[7440-31-5]	8,84	%		
	Zinc (Zn)	[7440-66-6]	6,3	%		