

**Copper base, Type: Low Alloy, RM, Chips 50 g**

Art. ID DH-SL02-02  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	0,003	%		
	Copper (Cu)	[7440-50-8]	98,15	%		
	Iron (Fe)	[7439-89-6]	0,911	%		
	Manganese (Mn)	[7439-96-5]	0,007	%		
	Nickel (Ni)	[7440-02-0]	0,034	%		
	Lead (Pb)	[7439-92-1]	0,139	%		
	Sulfur (S)	[7704-34-9]	0,037	%		
	Antimony (Sb)	[7440-36-0]	0,008	%		
	Tin (Sn)	[7440-31-5]	0,381	%		
	Zinc (Zn)	[7440-66-6]	0,229	%		