

**Copper base, Type: Low Alloy, RM, Chips 50 g**

Art. ID                      DH-SL02-04  
Unit                         50 g (chips)  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	12,51	%		
	Chromium (Cr)	[7440-47-3]	0,009	%		
	Copper (Cu)	[7440-50-8]	77	%		
	Iron (Fe)	[7439-89-6]	5,7	%		
	Manganese (Mn)	[7439-96-5]	0,057	%		
	Phosphorus (P)	[7723-14-0]	0,007	%		
	Lead (Pb)	[7439-92-1]	0,58	%		
	Antimony (Sb)	[7440-36-0]	0,336	%		
	Silicon (Si)	[7440-21-3]	0,22	%		
	Tin (Sn)	[7440-31-5]	2,16	%		
	Zinc (Zn)	[7440-66-6]	1,36	%		