

**Copper base, Type: Bronze, RM, Chips 50 g**

Art. ID                      DH-SL02-06  
Unit                         50 g (chips)  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,059	%		
	Arsenic (As)	[7440-38-2]	0,025	%		
	Copper (Cu)	[7440-50-8]	83,05	%		
	Iron (Fe)	[7439-89-6]	1,79	%		
	Manganese (Mn)	[7439-96-5]	0,044	%		
	Nickel (Ni)	[7440-02-0]	0,221	%		
	Phosphorus (P)	[7723-14-0]	0,017	%		
	Lead (Pb)	[7439-92-1]	0,891	%		
	Sulfur (S)	[7704-34-9]	0,059	%		
	Antimony (Sb)	[7440-36-0]	0,06	%		
	Silicon (Si)	[7440-21-3]	0,043	%		
	Tin (Sn)	[7440-31-5]	2,78	%		
	Zinc (Zn)	[7440-66-6]	10,89	%		