

Copper base, Type: Low Alloy, RM, Chips 50 g

Art. ID DH-SL02-08
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	4,15	%		
	Chromium (Cr)	[7440-47-3]	0,009	%		
	Copper (Cu)	[7440-50-8]	81,67	%		
	Iron (Fe)	[7439-89-6]	2,54	%		
	Manganese (Mn)	[7439-96-5]	0,711	%		
	Nickel (Ni)	[7440-02-0]	2,82	%		
	Phosphorus (P)	[7723-14-0]	0,027	%		
	Lead (Pb)	[7439-92-1]	1,31	%		
	Antimony (Sb)	[7440-36-0]	0,083	%		
	Silicon (Si)	[7440-21-3]	0,052	%		
	Tin (Sn)	[7440-31-5]	4,78	%		
	Zinc (Zn)	[7440-66-6]	1,85	%		