

Copper base, Type: Low Alloy, RM, Chips 50 g

Art. ID DH-SL02-09
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	87,07	%		
	Nickel (Ni)	[7440-02-0]	0,265	%		
	Lead (Pb)	[7439-92-1]	0,542	%		
	Tin (Sn)	[7440-31-5]	11,92	%		