

Iron Base, Steel - Low Alloy, RM, Chips 50 g

Art. ID DH-SE23-03
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Insoluble Al		0,05	%		
	Arsenic (As)	[7440-38-2]	0,003	%		
	Carbon (C)	[7440-44-0]	0,1671	%		
	Manganese (Mn)	[7439-96-5]	1,744	%		
	Phosphorus (P)	[7723-14-0]	0,0151	%		
	Sulfur (S)	[7704-34-9]	0,0112	%		
	Silicon (Si)	[7440-21-3]	0,577	%		