

Iron Base, Steel - Low Alloy, RM, Chips 50 g

Art. ID DH-SE23-29
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0436	%		
	Arsenic (As)	[7440-38-2]	0,002	%		
	Boron (B)	[7440-42-8]	<0,00009	%		
	Bismuth (Bi)	[7440-69-9]	<0,00002	%		
	Carbon (C)	[7440-44-0]	0,168	%		
	Calcium (Ca)	[7440-70-2]	0,0027	%		
	Cobalt (Co)	[7440-48-4]	0,005	%		
	Chromium (Cr)	[7440-47-3]	1,7057	%		
	Copper (Cu)	[7440-50-8]	0,0216	%		
	Manganese (Mn)	[7439-96-5]	0,3219	%		
	Molybdenum (Mo)	[7439-98-7]	0,4126	%		
	Nitrogen (N)	[7727-37-9]	0,0045	%		
	Niobium (Nb)	[7440-03-1]	0,00097	%		
	Nickel (Ni)	[7440-02-0]	2,3865	%		
	Phosphorus (P)	[7723-14-0]	0,0065	%		
	Lead (Pb)	[7439-92-1]	<0,00007	%		
	Sulfur (S)	[7704-34-9]	0,00067	%		
	Antimony (Sb)	[7440-36-0]	0,00053	%		
	Silicon (Si)	[7440-21-3]	0,257	%		
	Tin (Sn)	[7440-31-5]	0,00052	%		
	Titanium (Ti)	[7440-32-6]	0,0027	%		
	Vanadium (V)	[7440-62-2]	0,0025	%		