

Iron Base, Steel - High Alloy, RM, Chips 50 g

Art. ID DH-SE24-01
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Soluble Al		0,1013	%		
	Arsenic (As)	[7440-38-2]	0,0038	%		
	Carbon (C)	[7440-44-0]	0,0048	%		
	Calcium (Ca)	[7440-70-2]	0,0002	%		
	Cobalt (Co)	[7440-48-4]	11,46	%		
	Chromium (Cr)	[7440-47-3]	0,0405	%		
	Copper (Cu)	[7440-50-8]	0,0161	%		
	Manganese (Mn)	[7439-96-5]	0,0203	%		
	Molybdenum (Mo)	[7439-98-7]	4,64	%		
	Niobium (Nb)	[7440-03-1]	0,0023	%		
	Nickel (Ni)	[7440-02-0]	18,85	%		
	Phosphorus (P)	[7723-14-0]	0,0013	%		
	Lead (Pb)	[7439-92-1]	<0,0001	%		
	Sulfur (S)	[7704-34-9]	0,0026	%		
	Antimony (Sb)	[7440-36-0]	0,0003	%		
	Silicon (Si)	[7440-21-3]	0,0222	%		
	Tin (Sn)	[7440-31-5]	0,0016	%		
	Titanium (Ti)	[7440-32-6]	1,28	%		
	Vanadium (V)	[7440-62-2]	0,0029	%		