

## **Proficiency Testing - PCBs in Oil**

Art. ID                      ERA-817  
Unit                         each

### Description

One 10 mL flame-sealed ampule ready to analyze. Use with EPA method 8082, or other applicable method. Contains a different Aroclor, randomly selected at 10-50 mg/kg.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aroclor 1016	[12674-11-2]	10 - 50	mg/kg		
	Aroclor 1221	[11104-28-2]	12 - 50	mg/kg		
	Aroclor 1232	[11141-16-5]	12 - 50	mg/kg		
	Aroclor 1242	[53469-21-9]	10 - 50	mg/kg		
	Aroclor 1248	[12672-29-6]	12 - 50	mg/kg		
	Aroclor 1254	[11097-69-1]	10 - 50	mg/kg		
	Aroclor 1260	[11096-82-5]	10 - 50	mg/kg		