

Assurance Grade Arsenic (+3) Speciation Standard for ICP and ICP-MS - Arsenic-III (As+3) 1000 mg/L in 2% HCl

Art. ID LM24-SP-SPEC-AS3M
Unit 30 mL
Deliverydetails No Dangerous Good UN1789 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|------|--------|--------|
| | Arsenic-III (As+3) | | 1000 | mg/L | | |