

Assurance Grade Selenium (+4) Speciation Standard for ICP and ICP-MS - Selenium-IV (Se+4) 1000 mg/L in 2% HNO₃

Art. ID LM24-SP-SPEC-SE4M
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Selenium-IV (Se+4)		1000	mg/L		