

Gonyautoxin-2 and3 (GTX2and3) at a concentration of 146,1 µmol/L in 0,5 mM hydrochloric acid

Art. ID NRC-CRM-GTX2and3-d
Unit 0,5 mL/ampoule
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
GTX2	Gonyautoxin 2	[60508-89-6]	102,6 ± 14			
GTX3	Gonyautoxin 3	[60537-65-7]	43,5 ± 2,9			
GTX2 + GTX3	Gonyautoxin 2 (GTX2) + Gonyautoxin 3 (GTX3)	[60508-89-6]/[6 0537-65-7]	146,1 ± 6,9			