

Corrosion-Resistant Cast Iron, Chill Cast, Certified Reference Material (CRM), disc 40 mm x 13 mm

Art. ID MBH-11X0331.3H
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,08	%		
	Silicon (Si)	[7440-21-3]	1,83	%		
	Sulfur (S)	[7704-34-9]	0,056	%		
	Phosphorus (P)	[7723-14-0]	0,0459	%		
	Manganese (Mn)	[7439-96-5]	0,756	%		
	Nickel (Ni)	[7440-02-0]	16,51	%		
	Chromium (Cr)	[7440-47-3]	1,53	%		
	Molybdenum (Mo)	[7439-98-7]	0,0381	%		
	Copper (Cu)	[7440-50-8]	7,01	%		
	Cobalt (Co)	[7440-48-4]	0,126	%		
	Vanadium (V)	[7440-62-2]	0,0386	%		
	Niobium (Nb)	[7440-03-1]	0,071	%		
	Titanium (Ti)	[7440-32-6]	0,075	%		
	Aluminium (Al)	[7429-90-5]	0,0323	%		
	Tin (Sn)	[7440-31-5]	0,0098	%		
	Lead (Pb)	[7439-92-1]	0,0088	%		