

ICP-MS Instrument Calibration Standard 2 in 5% HNO₃ / tr Tart.

Art. ID LM24-SP-CL-CAL-2
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]		100	µg/mL	
	Antimony (Sb)	[7440-36-0]		100	µg/mL	
	Arsenic (As)	[7440-38-2]		100	µg/mL	
	Barium (Ba)	[7440-39-3]		100	µg/mL	
	Beryllium (Be)	[7440-41-7]		100	µg/mL	
	Cadmium (Cd)	[7440-43-9]		100	µg/mL	
	Calcium (Ca)	[7440-70-2]		100	µg/mL	
	Chromium (Cr)	[7440-47-3]		100	µg/mL	
	Cobalt (Co)	[7440-48-4]		100	µg/mL	
	Copper (Cu)	[7440-50-8]		100	µg/mL	
	Iron (Fe)	[7439-89-6]		100	µg/mL	
	Lead (Pb)	[7439-92-1]		100	µg/mL	
	Magnesium (Mg)	[7439-95-4]		100	µg/mL	
	Manganese (Mn)	[7439-96-5]		100	µg/mL	
	Molybdenum (Mo)	[7439-98-7]		100	µg/mL	
	Nickel (Ni)	[7440-02-0]		100	µg/mL	
	Potassium (K)	[7440-09-7]		100	µg/mL	
	Selenium (Se)	[7782-49-2]		100	µg/mL	
	Silver (Ag)	[7440-22-4]		100	µg/mL	
	Sodium (Na)	[7440-23-5]		100	µg/mL	
	Strontium (Sr)	[7440-24-6]		100	µg/mL	
	Thallium (Tl)	[7440-28-0]		100	µg/mL	
	Tin (Sn)	[7440-31-5]		100	µg/mL	
	Titanium (Ti)	[7440-32-6]		100	µg/mL	
	Vanadium (V)	[7440-62-2]		100	µg/mL	
	Zinc (Zn)	[7440-66-6]		100	µg/mL	