

Genetically modified DAS-81910-7 Cotton (0,1%)

Art. ID ERM-BF440c
Unit 1 g cotton powder
Deliverydetails No Dangerous Good /not restricted

Description

ERM-BF440c is intended to be used for calibration or quality control of methods for the identification and quantification of genetically modified DAS-81910-7 cotton in food and feed.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DAS-81910-7 cotton		1 ± 0,08	g/kg		