

Genetically modified DAS-81910-7 Cotton (1%)

Art. ID ERM-BF440d
Unit 1 g cotton powder
Deliverydetails No Dangerous Good /not restricted

Description

ERM-BF440d is intended to be used for calibration or quality control of methods for the identification and quantification of genetically modified DAS-81910-7 cotton in food and feed.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|------|--------|--------|
| | DAS-81910-7 cotton | | 10 ± 0,8 | g/kg | | |