

Custom Solution - 5 components 100 mg/L each of in HCl 5%

Art. ID LM24-CUS-70224
Unit 2 x 25 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 II E2

Description

Stored under normal laboratory conditions. Conservation 24 months. (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iridium (Ir)	[7439-88-5]	100	mg/L		
	Palladium (Pd)	[7440-05-3]	100	mg/L		
	Platinum (Pt)	[7440-06-4]	100	mg/L		
	Rhodium (Rh)	[7440-16-6]	100	mg/L		
	Ruthenium (Ru)	[7440-18-8]	100	mg/L		