

Calibration Blank (Anatel PAT700)

Art. ID REA-ISTOC1001
Unit 60 mL
Deliverydetails No Dangerous Good /not restricted

Description

Reagecon's Calibration Blank Total Organic Carbon (TOC) Standard is a product of extremely high specification and purity. The product is produced from water that is prepared by a special proprietary process. The product is presented in bottles that are manufactured, cleaned and leached, by a special proprietary process specifically for the product. Current Certificates are 17025 approved / Traceable to a NIST standard

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------|------------|---------------------|------|--------|--------|
| | Total Organic Carbon (TOC) | | Blank | | | |