

## **USP Reagent Water System Suitability Standard (Anatel PAT700)**

Art. ID REA-ISTOC1005  
Unit 60 mL  
Deliverydetails No Dangerous Good /not restricted

### **Description**

Reagecon's USP Reagent Water Total Organic Carbon (TOC) Standard is a product of extremely high specification and purity. The product is produced from water that is prepared by a special proprietary process. The product is presented in bottles that are manufactured, cleaned and leached, by a special proprietary process specifically for the product. Current Certificates are 17025 approved / Traceable to a NIST standard

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Organic Carbon (TOC)					