

Calibration Kit, 2 x blank, 2 x 0.1 mg/L, 2 x 0.25 mg/L , 2 x 0.5 mg/L , 0.75 mg/L, 1 mg/L C NIST KHP (Shimadzu)

| | |
|------------------|-----------------------------------|
| Art. ID | REA-ISTOC1044 |
| Unit | Kit (10 x 40 mL) |
| Delivery details | No Dangerous Good /not restricted |

Description

The item has a very short shelf life / Current Certificates are 17025 approved / Traceable to a NIST standard