

### **Anatel A643 Calibration Blank (Anatel A643 & Anatel TOC600)**

Art. ID	REA-ISTOC1079
Unit	60 mL
Deliverydetails	No Dangerous Good /not restricted

#### Description

Reagecon's Calibration Blank Total Organic Carbon (TOC) Standard is a product of extremely high specification and purity. The product is produced from water that is prepared by a special proprietary process. The product is presented in bottles that are manufactured, cleaned and leached, by a special proprietary process specifically for the product.. Current Certificates are 17025 approved / Traceable to a NIST standard

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Organic Carbon (TOC)		Blank			