

**\*\*\*not available anymore\*\*\*Phosphonic acid-18O3 - approx. 5% solution in water-18O with an isotope enrichment (\*) of 97%**

Art. ID LM24-AC-20250  
Unit on request  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2834 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphonic acid-18O3	[13598-36-2] (unlabelled)				