

**\*\*\*not available anymore\*\*\*Phosphonic acid-18O3 - approx. 5% solution in water-18O with an isotope enrichment (\*) of 97%**

Art. ID LM24-AC-20250  
Unit on request  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2834 III E1

| Text/Information | Analyte/Parameter    | CAS number                | Concentration/Value | Unit | Method | Source |
|------------------|----------------------|---------------------------|---------------------|------|--------|--------|
|                  | Phosphonic acid-18O3 | [13598-36-2] (unlabelled) |                     |      |        |        |