

Custom Solution - 8 components in 7% HNO₃

Art. ID LM24-CUS-70232
Unit 2 x 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Sulfur (S)	[7704-34-9]	2000	mg/L		
	Potassium (K)	[7440-09-7]	500	mg/L		
	Silicon (Si)	[7440-21-3]	250	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Calcium (Ca)	[7440-70-2]	500	mg/L		
	Sodium (Na)	[7440-23-5]	1000	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		