

## **Custom Solution - 6 components in water**

Art. ID LM24-CUS-70234  
Unit 25 mL  
Deliverydetails No Dangerous Good /not restricted

### Description

Stored under normal laboratory conditions. Conservation 12 months.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fluoride (F-)		400	mg/L		
	Nitrite (NO <sub>2</sub> -)		250	mg/L		
	Nitrate (NO <sub>3</sub> -)		400	mg/L		
	Phosphate (PO <sub>4</sub> <sup>3-</sup> )		250	mg/L		
	Chloride (Cl-)		4000	mg/L		
	Sulfate (SO <sub>4</sub> <sup>2-</sup> )		20000	mg/L		