

Custom Solution - 2 components in water

Art. ID LM24-CUS-70235
Unit 25 mL
Deliverydetails No Dangerous Good /not restricted

Description

Stored under normal laboratory conditions. Conservation 12 months.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|------|--------|--------|
| | Phosphates (PO43-) | | 1338 | mg/L | | |
| | Potassium (K+) | | 830 | mg/L | | |