

Potassium ICP standard traceable to SRM from NIST KNO₃ in HNO₃ 2-3% 10000 mg/L K Certipur®

Art. ID SAF-1703770100
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Potassium (K) | [7440-09-7] | 10000 | mg/L | | |