

### **Custom anion standard - Metals in water**

Art. ID LM24-CUS-70013  
Unit 2 x 25mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	20	µg/L		
	Boron (B)	[7440-42-8]	100	µg/L		
	Cadmium (Cd)	[7440-43-9]	0,5	µg/L		
	Chromium (Cr)	[7440-47-3]	5	µg/L		
	Copper (Cu)	[7440-50-8]	200	µg/L		
	Iron (Fe)	[7439-89-6]	20	µg/L		
	Manganese (Mn)	[7439-96-5]	5	µg/L		
	Nickel (Ni)	[7440-02-0]	2	µg/L		
	Lead (Pb)	[7439-92-1]	1	µg/L		
	Antimony (Sb)	[7440-36-0]	1,25	µg/L		
	Arsenic (As)	[7440-38-2]	1	µg/L		
	Selenium (Se)	[7782-49-2]	1	µg/L		