

**Austenitic (Ni-RESIST), Cast Iron, RM, Block 48 mm x 42 mm x 12 mm**

Art. ID                      BAS-NIRM 5/1  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,95	%		
	Cerium (Ce)	[7440-45-1]	0,016	%		
	Chromium (Cr)	[7440-47-3]	0,51	%		
	Copper (Cu)	[7440-50-8]	0,21	%		
	Magnesium (Mg)	[7439-95-4]	0,055	%		
	Manganese (Mn)	[7439-96-5]	1,01	%		
	Niobium (Nb)	[7440-03-1]	0,15	%		
	Nickel (Ni)	[7440-02-0]	21,7	%		
	Phosphorus (P)	[7723-14-0]	0,103	%		
	Sulfur (S)	[7704-34-9]	0,005	%		
	Silicon (Si)	[7440-21-3]	1,50	%		