

High Manganese Steel, Iron Base (Cast), CRM, Block 48 mm x 42 mm x 12 mm

Art. ID BAS-SS-CRM 493/3
Unit Disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,035	%		
	Carbon (C)	[7440-44-0]	0,819	%		
	Chromium (Cr)	[7440-47-3]	0,259	%		
	Copper (Cu)	[7440-50-8]	0,017	%		
	Manganese (Mn)	[7439-96-5]	11,15	%		
	Molybdenum (Mo)	[7439-98-7]	1,04	%		
	Nitrogen (N)	[7727-37-9]	0,024	%		
	Nickel (Ni)	[7440-02-0]	3,24	%		
	Phosphorus (P)	[7723-14-0]	0,12	%		
	Sulfur (S)	[7704-34-9]	0,009	%		
	Silicon (Si)	[7440-21-3]	0,861	%		
	Vanadium (V)	[7440-62-2]	0,025	%		