

EPA Method 624 Purgeable Mixture in Methanol

Art. ID LM24-PMX-110-1
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	100	µg/mL		
	Bromodichloromethane	[75-27-4]	100	µg/mL		
	Bromoform (Tribromomethane)	[75-25-2]	100	µg/mL		
	Bromomethane	[74-83-9]	100	µg/mL		
	Carbon tetrachloride (Tetrachloromethane)	[56-23-5]	100	µg/mL		
	Chlorobenzene	[108-90-7]	100	µg/mL		
	Chloroethane	[75-00-3]	100	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	100	µg/mL		
	Chloromethane	[74-87-3]	100	µg/mL		
	Dibromochloromethane	[124-48-1]	100	µg/mL		
	1,2-Dichlorobenzene	[95-50-1]	100	µg/mL		
	1,3-Dichlorobenzene	[541-73-1]	100	µg/mL		
	1,4-Dichlorobenzene	[106-46-7]	100	µg/mL		
	1,1-Dichloroethane	[75-34-3]	100	µg/mL		
	1,2-Dichloroethane	[107-06-2]	100	µg/mL		
	1,1-Dichloroethene	[75-35-4]	100	µg/mL		
	trans-1,2-Dichloroethene	[156-60-5]	100	µg/mL		
	1,2-Dichloropropane	[78-87-5]	100	µg/mL		
	cis-1,3-Dichloropropene	[10061-01-5]	100	µg/mL		
	trans-1,3-Dichloropropene	[10061-02-6]	100	µg/mL		
	Ethylbenzene	[100-41-4]	100	µg/mL		
	Methylene chloride (Dichloromethane)	[75-09-2]	100	µg/mL		
	1,1,2,2-Tetrachloroethane	[79-34-5]	100	µg/mL		
	Tetrachloroethene	[127-18-4]	100	µg/mL		

Toluene	[108-88-3]	100	µg/mL
1,1,1-Trichloroethane	[71-55-6]	100	µg/mL
1,1,2-Trichloroethane	[79-00-5]	100	µg/mL
Trichloroethene	[79-01-6]	100	µg/mL
Trichlorofluoromethane	[75-69-4]	100	µg/mL
Vinyl chloride	[75-01-4]	100	µg/mL