

Iron Base, Cast Iron (vermicular iron, CGI), CRM, Disc 40 mm x 18 mm

Art. ID CZ-02033-3C
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,024 | % | | |
| | Boron (B) | [7440-42-8] | 0,0044 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,002 | % | | |
| | Carbon (C) | [7440-44-0] | 3,68 | % | | |
| | Cerium (Ce) | [7440-45-1] | 0,013 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,026 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,1 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,421 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,006 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,333 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,49 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,04 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,026 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,005 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,007 | % | | |
| | Silicon (Si) | [7440-21-3] | 2,15 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,009 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,021 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,016 | % | | |
| | Tungsten (W) | [7440-33-7] | 0,003 | % | | |