

Custom Solution - 6 elements in nitric acid (5%)

Art. ID LM24-CUS-70020
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Expiration 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	300	mg/L		
	Molybdenum (Mo)	[7439-98-7]	300	mg/L		
	Silicon (Si)	[7440-21-3]	300	mg/L		
	Phosphorus (P)	[7723-14-0]	2500	mg/L		
	Sulfur (S)	[7704-34-9]	17000	mg/L		