

Al-0.1%Au - Aluminium gold alloy certified for the mass fraction of gold (IRMM)

Art. ID ERM-EB530B
Unit wire
Deliverydetails No Dangerous Good /not restricted

Description

Certified value. The ERM-EB530B unit is a 1 meter lengtht wire with diameter of 0,500 mm (relative tolerance $\pm 5\%$).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
1005 $\pm$ 7 mg/kg	Aluminium (Al)	[7429-90-5]	0,1005	%		