

**Cast Iron, CRM, Truncated cone ~35 mm Ø x ~20 mm**

Art. ID VS-ChG31  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	0,068	%		
	Carbon (C)	[7440-44-0]	3,19	%		
	Chromium (Cr)	[7440-47-3]	0,156	%		
	Copper (Cu)	[7440-50-8]	0,281	%		
	Manganese (Mn)	[7439-96-5]	0,97	%		
	Molybdenum (Mo)	[7439-98-7]	0,0069	%		
	Phosphorus (P)	[7723-14-0]	0,047	%		
	Sulfur (S)	[7704-34-9]	0,043	%		
	Silicon (Si)	[7440-21-3]	1,6	%		
	Tin (Sn)	[7440-31-5]	0,013	%		
	Titanium (Ti)	[7440-32-6]	0,0063	%		
	Vanadium (V)	[7440-62-2]	0,0035	%		