

**Cast Iron - High Alloy, CRM, Truncated cone ~35 mm Ø x ~17 mm**

Art. ID VS-ChG43  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,6	%		
	Chromium (Cr)	[7440-47-3]	22,79	%		
	Copper (Cu)	[7440-50-8]	0,24	%		
	Manganese (Mn)	[7439-96-5]	1,26	%		
	Molybdenum (Mo)	[7439-98-7]	~0,02	%		
	Nickel (Ni)	[7440-02-0]	0,28	%		
	Phosphorus (P)	[7723-14-0]	0,052	%		
	Sulfur (S)	[7704-34-9]	0,015	%		
	Silicon (Si)	[7440-21-3]	3,87	%		
	Titanium (Ti)	[7440-32-6]	0,036	%		
	Vanadium (V)	[7440-62-2]	0,028	%		