

**Cast Iron, CRM, Pig iron, Truncated cone ~38 mm Ø x ~38 mm**

Art. ID VS-ChG5/9  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	~0,003	%		
	Carbon (C)	[7440-44-0]	3,51	%		
	Chromium (Cr)	[7440-47-3]	0,307	%		
	Copper (Cu)	[7440-50-8]	0,037	%		
	Manganese (Mn)	[7439-96-5]	0,6	%		
	Phosphorus (P)	[7723-14-0]	0,104	%		
	Sulfur (S)	[7704-34-9]	0,036	%		
	Silicon (Si)	[7440-21-3]	0,84	%		
	Titanium (Ti)	[7440-32-6]	~0,1	%		
	Vanadium (V)	[7440-62-2]	0,441	%		