

**Tin White Metal (Pewter Babbitt Alloy), Cast, Certified Reference Material (CRM), disc 40 mm x 15 mm**

Art. ID MBH-73XSC2B  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0318	%		
	Arsenic (As)	[7440-38-2]	0,0355	%		
	Bismuth (Bi)	[7440-69-9]	0,602	%		
	Cadmium (Cd)	[7440-43-9]	0,0282	%		
	Cobalt (Co)	[7440-48-4]	0,0217	%		
	Copper (Cu)	[7440-50-8]	1,028	%		
	Iron (Fe)	[7439-89-6]	0,0037	%		
	Indium (In)	[7440-74-6]	0,0099	%		
	Nickel (Ni)	[7440-02-0]	0,0156	%		
	Lead (Pb)	[7439-92-1]	0,0457	%		
	Antimony (Sb)	[7440-36-0]	3,05	%		
	Zinc (Zn)	[7440-66-6]	~0,0007	%		