

## **Wastewater - Suspended solids (SS)**

Art. ID LM-VKI-QC-WW6-50  
Unit 50 ampoules  
Deliverydetails No Dangerous Good /not restricted

### **Description**

Produced and certified under DANAK accreditation no. 535 according to ISO 17034 /// LM-VKI-QC-WW6 is a bottle containing 3.00 g concentrate from which 83 mL reference sample is prepared /// WW6 consists of a bottle of reference material for suspended solids with a concentration of 236 mg/L, but may also be diluted to 23,6 or 2.36 mg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Suspended solids (SS)		236	mg/L		