

**Wastewater/Fresh recipient water/Drinking water and Groundwater -
Al-As-Cr-Cu-Mo-Ni-Sn-V-Zn: ppb-level**

Art. ID LM-VKI-QC-METAL-LL1
Unit 1 ampoule
Deliverydetails No Dangerous Good

Description

1L prepared reference material per ampoule

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	200	µg/L		
	Arsenic (As)	[7440-38-2]	30	µg/L		
	Chromium (Cr)	[7440-47-3]	20	µg/L		
	Copper (Cu)	[7440-50-8]	20	µg/L		
	Molybdenum (Mo)	[7439-98-7]	20	µg/L		
	Nickel (Ni)	[7440-02-0]	20	µg/L		
	Tin (Sn)	[7440-31-5]	20	µg/L		
	Vanadium (V)	[7440-62-2]	20	µg/L		
	Zinc (Zn)	[7440-66-6]	50	µg/L		