

**Wastewater/Fresh recipient water/Drinking water and Groundwater -  
Ag-Ba-Cd-Co-Fe-Mn-Pb-Sb-Se-Sr: ppb-level**

Art. ID LM-VKI-QC-METAL-LL2  
Unit 1 ampoule  
Deliverydetails No Dangerous Good UN2031 II E2

Description

1L prepared reference material per ampoule

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	5	µg/L		
	Barium (Ba)	[7440-39-3]	100	µg/L		
	Cadmium (Cd)	[7440-43-9]	2	µg/L		
	Cobalt (Co)	[7440-48-4]	50	µg/L		
	Iron (Fe)	[7439-89-6]	200	µg/L		
	Manganese (Mn)	[7439-96-5]	50	µg/L		
	Lead (Pb)	[7439-92-1]	20	µg/L		
	Antimony (Sb)	[7440-36-0]	50	µg/L		
	Selenium (Se)	[7782-49-2]	100	µg/L		
	Strontium (Sr)	[7440-24-6]	50	µg/L		